

Buzz - Upholstery

Care 6

MAINTENANCE
INSTRUCTION FOR
FABRICS WITH
CRYPTON® FINISH

Cleaning Code W: Clean only with water-based shampoo or foam upholstery cleaner. Do not over wet. Do not use dry cleaning solvents to spot clean. Pile fabrics may require brushing to restore appearance. Cushion covers should not be removed and laundered. Diluted bleach solution can be used, if needed for disinfecting.

- Most spills on Crypton can easily be cleaned with a towel. Just blot up the spill with a clean cloth.
- Dried stains or more difficult stains can be easily cleaned with a soap-and-water mixture.
- Simply combine 1 part of an enzyme powder detergent such as Tide® or Cheer® with 5 parts water. Apply the solution to the stain and wait a minute or two to allow the detergent to begin removing the stain. Then lightly rub/scrub the area and rinse thoroughly with clean water.
- After cleaning, always rinse with fresh water to remove any detergent or cleaning agent residue from the surface of the fabric.
- For Dried Heavy Liquids: Slide a dull, flat edge under the staining agent, hold a clean cloth in the opposite hand and 'lift' the staining agent from the fabric onto the clean cloth. If 'heavy liquid' has completely dried, try 'peeling' the stain off with your fingers. If stain does not 'peel' off completely, continue to try to remove as much of the staining agent as possible soap and water.
- If stain is oil-based, pre-treat the area with rubbing alcohol. Liberally apply the alcohol and rub into stain. Allow alcohol to sit on the stain for 1 to 2 minutes, blot with a clean cloth. Be certain to rinse thoroughly with clean water and dry with a clean cloth.
- For disinfecting, a 10% bleach solution may be used on this fabric.